



FCER1A Antibody, FITC conjugated

Product Code	CSB-PA008532LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12319
Immunogen	Recombinant Human High affinity immunoglobulin epsilon receptor subunit alpha protein (26-205AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	High affinity immunoglobulin epsilon receptor subunit alpha (Fc-epsilon RI-alpha) (FcERI) (IgE Fc receptor subunit alpha), FCER1A, FCE1A
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCER1A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.