





FOXG1 Antibody

Product Code	CSB-PA008815HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55316
Immunogen	Recombinant Human Forkhead box protein G1 protein (183-292AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein G1 (Brain factor 1) (BF-1) (BF1) (Brain factor 2) (BF-2) (BF2) (hBF-2) (Forkhead box protein G1A) (Forkhead box protein G1B) (Forkhead box protein G1C) (Forkhead-related protein FKHL1) (HFK1) (Forkhead-related protein FKHL2) (HFK2) (Forkhead-related protein FKHL3) (HFK3), FOXG1, FKHL1 FKHL2 FKHL3 FKHL4 FOXG1A FOXG1B FOXG1C
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FOXG1

Image



Immunohistochemistry of paraffin-embedded human heart tissue using CSB-PA008815HA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA008815HA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)