



MTOR Antibody, HRP conjugated

Product Code	CSB-PA008968HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P42345
Immunogen	Recombinant Human Serine/threonine-protein kinase mTOR protein (1223-1290AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase mTOR (EC 2.7.11.1) (FK506-binding protein 12-rapamycin complex-associated protein 1) (FKBP12-rapamycin complex-associated protein) (Mammalian target of rapamycin) (mTOR) (Mechanistic target of rapamycin) (Rapamycin and FKBP12 target 1) (Rapamycin target protein 1), MTOR, FRAP FRAP1 FRAP2 RAFT1 RAPT1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MTOR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.