



G6PD Antibody, FITC conjugated

Product Code	CSB-PA009121LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P11413
Immunogen	Recombinant Human Glucose-6-phosphate 1-dehydrogenase protein (80-185AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glucose-6-phosphate 1-dehydrogenase (G6PD) (EC 1.1.1.49), G6PD
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	G6PD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.