



GAA Antibody, Biotin conjugated

Product Code	CSB-PA009125LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10253
Immunogen	Recombinant Human Lysosomal alpha-glucosidase protein (601-952AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysosomal alpha-glucosidase (EC 3.2.1.20) (Acid maltase) (Aglucosidase alfa) [Cleaved into: 76 kDa lysosomal alpha-glucosidase; 70 kDa lysosomal alpha-glucosidase], GAA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GAA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.