



GLUL Antibody, HRP conjugated

Product Code	CSB-PA009553LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P15104
Immunogen	Recombinant Human Glutamine synthetase protein (86-373AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutamine synthetase (GS) (EC 6.3.1.2) (Glutamate decarboxylase) (EC 4.1.1.15) (Glutamate--ammonia ligase), GLUL, GLNS
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GLUL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.