





GRID2IP Antibody

Product Code	CSB-PA009904LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A4D2P6
Immunogen	Recombinant Human Delphilin protein (35-156AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Delphilin (Glutamate receptor, ionotropic, delta 2-interacting protein 1), GRID2IP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GRID2IP


Image



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA009904LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA009904LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Immunofluorescent analysis of MCF-7 cells using CSB-PA009904LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)