



GUCY1A1 Antibody, Biotin conjugated

Product Code	CSB-PA010050LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02108
Immunogen	Recombinant Human Guanylate cyclase soluble subunit alpha-1 protein (20-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Guanylate cyclase soluble subunit alpha-1 (GCS-alpha-1) (EC 4.6.1.2) (Guanylate cyclase soluble subunit alpha-3) (GCS-alpha-3) (Soluble guanylate cyclase large subunit), GUCY1A1, GUC1A3 GUCSA3 GUCY1A3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GUCY1A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.