



HFM1 Antibody, Biotin conjugated

Product Code	CSB-PA010325LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A2PYH4
Immunogen	Recombinant Human Probable ATP-dependent DNA helicase HFM1 protein (121-420AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable ATP-dependent DNA helicase HFM1 (EC 3.6.4.12) (SEC63 domain-containing protein 1), HFM1, SEC3D1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	HFM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.