



HOXD9 Antibody, Biotin conjugated

Product Code	CSB-PA010689LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28356
Immunogen	Recombinant Human Homeobox protein Hox-D9 protein (38-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Hox-D9 (Homeobox protein Hox-4C) (Homeobox protein Hox-5.2), HOXD9, HOX4C
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HOXD9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.