



HS6ST1 Antibody, HRP conjugated

Product Code	CSB-PA010756LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O60243
Immunogen	Recombinant Human Heparan-sulfate 6-O-sulfotransferase 1 protein (3-133AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heparan-sulfate 6-O-sulfotransferase 1 (HS6ST-1) (EC 2.8.2.-), HS6ST1, HS6ST
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HS6ST1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.