



Hspa5 Antibody, HRP conjugated

Product Code	CSB-PA010827YB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20029
Immunogen	Recombinant Mouse Endoplasmic reticulum chaperone BiP protein (20-655AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoplasmic reticulum chaperone BiP (EC 3.6.4.10) (78 kDa glucose-regulated protein) (GRP-78) (Binding-immunoglobulin protein) (BiP) (Heat shock protein 70 family protein 5) (HSP70 family protein 5) (Heat shock protein family A member 5) (Immunoglobulin heavy chain-binding protein), Hspa5, Grp78
Immunogen Species	Mus musculus (Mouse)
Research Area	Tags & Cell Markers
Target Names	Hspa5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.