



# Hspa5 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA010827YC01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P20029
<b>Immunogen</b>	Recombinant Mouse Endoplasmic reticulum chaperone BiP protein (20-655AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Endoplasmic reticulum chaperone BiP (EC 3.6.4.10) (78 kDa glucose-regulated protein) (GRP-78) (Binding-immunoglobulin protein) (BiP) (Heat shock protein 70 family protein 5) (HSP70 family protein 5) (Heat shock protein family A member 5) (Immunoglobulin heavy chain-binding protein), Hspa5, Grp78
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Tags & Cell Markers
<b>Target Names</b>	Hspa5
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.