





CHD4 Antibody

Product Code	CSB-PA01097A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14839
Immunogen	Recombinant Human Chromodomain-helicase-DNA-binding protein 4 protein (1-219AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chromodomain-helicase-DNA-binding protein 4 (CHD-4) (EC 3.6.4.12) (ATP-dependent helicase CHD4) (Mi-2 autoantigen 218 kDa protein) (Mi2-beta), CHD4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CHD4


Image



Immunohistochemistry of paraffin-embedded human liver tissue using CSB-PA01097A0Rb at dilution of 1:100



Immunohistochemistry of paraffin-embedded human melanoma using CSB-PA01097A0Rb at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA01097A0Rb at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.