



IFITM5 Antibody, Biotin conjugated

Product Code	CSB-PA011028HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NNB3
Immunogen	Recombinant Human Interferon-induced transmembrane protein 5 protein (1-36AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interferon-induced transmembrane protein 5 (Bone-restricted interferon-induced transmembrane protein-like protein) (BRIL) (Dispanin subfamily A member 1) (DSPA1), IFITM5
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	IFITM5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.