



IL13 Antibody, Biotin conjugated

Product Code	CSB-PA011590LD01HU
Abbreviation	Interleukin-13
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P35225
Immunogen	Recombinant Human Interleukin-13 protein (30-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	Cytokine (PubMed:8096327, PubMed:8097324). Inhibits inflammatory cytokine production (PubMed:8096327). Synergizes with IL2 in regulating interferon-gamma synthesis (PubMed:8096327). May be critical in regulating inflammatory and immune responses (PubMed:8096327, PubMed:8097324). Positively regulates IL31RA expression in macrophages (By similarity).
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-13 (IL-13), IL13, NC30
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	IL13