




INPP4B Antibody

Product Code	CSB-PA011732HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15327
Immunogen	Recombinant Human Type II inositol 3,4-bisphosphate 4-phosphatase protein (1-53AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Type II inositol 3,4-bisphosphate 4-phosphatase (EC 3.1.3.66) (Inositol polyphosphate 4-phosphatase type II), INPP4B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INPP4B

Image



Immunohistochemistry of paraffin-embedded human pancreatic cancer using CSB-PA011732HA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.