



INPP4B Antibody, Biotin conjugated

Product Code	CSB-PA011732HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15327
Immunogen	Recombinant Human Type II inositol 3,4-bisphosphate 4-phosphatase protein (1-53AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Type II inositol 3,4-bisphosphate 4-phosphatase (EC 3.1.3.66) (Inositol polyphosphate 4-phosphatase type II), INPP4B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INPP4B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.