



INTS9 Antibody, Biotin conjugated

Product Code	CSB-PA011767LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NV88
Immunogen	Recombinant Human Integrator complex subunit 9 protein (368-658AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrator complex subunit 9 (Int9) (Protein related to CPSF subunits of 74 kDa) (RC-74), INTS9, RC74
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	INTS9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.