



ACOD1 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA011825LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | A6NK06 |
| Immunogen | Recombinant Human Cis-aconitate decarboxylase protein (57-257AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Cis-aconitate decarboxylase (CAD) (EC 4.1.1.6) (Aconitate decarboxylase) (Aconitate decarboxylase 1) (Cis-aconitic acid decarboxylase) (Immune-responsive gene 1 protein), ACOD1, IRG1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | ACOD1 |