



ACOD1 Antibody, Biotin conjugated

Product Code	CSB-PA011825LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NK06
Immunogen	Recombinant Human Cis-aconitate decarboxylase protein (57-257AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cis-aconitate decarboxylase (CAD) (EC 4.1.1.6) (Aconitate decarboxylase) (Aconitate decarboxylase 1) (Cis-aconitic acid decarboxylase) (Immune-responsive gene 1 protein), ACOD1, IRG1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ACOD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.