



# ITGA4 Antibody, HRP conjugated

Product Code	CSB-PA011867LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13612
Immunogen	Recombinant Human Integrin alpha-4 protein (376-558AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-4 (CD49 antigen-like family member D) (Integrin alpha-IV) (VLA-4 subunit alpha) (CD antigen CD49d), ITGA4, CD49D
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ITGA4