



ITGA4 Antibody, Biotin conjugated

Product Code	CSB-PA011867LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13612
Immunogen	Recombinant Human Integrin alpha-4 protein (376-558AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Integrin alpha-4 (CD49 antigen-like family member D) (Integrin alpha-IV) (VLA-4 subunit alpha) (CD antigen CD49d), ITGA4, CD49D
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	ITGA4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.