



Lgals4 Antibody, HRP conjugated

Product Code	CSB-PA012889HB01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P38552
Immunogen	Recombinant Rat Galectin-4 protein (1-324AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Galectin-4 (Gal-4) (L-36 lactose-binding protein) (L36LBP) (Lactose-binding lectin 4), Lgals4
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Tags & Cell Markers
Target Names	Lgals4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.