




LMO2 Antibody

Product Code	CSB-PA013009HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25791
Immunogen	Recombinant Human Rhombotin-2 protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rhombotin-2 (Cysteine-rich protein TTG-2) (LIM domain only protein 2) (LMO-2) (T-cell translocation protein 2), LMO2, RBTN2 RBTN1 RHOM2 TTG2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LMO2

Image



Western blot


All lanes: LMO2 antibody at 8 µg/ml + Mouse spleen tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 19, 26, 6 kDa

Observed band size: 19 kDa



Immunohistochemistry of paraffin-embedded human tonsil tissue using CSB-PA013009HA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.