



LMO2 Antibody, HRP conjugated

Product Code	CSB-PA013009HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P25791
Immunogen	Recombinant Human Rhombotin-2 protein (1-100AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rhombotin-2 (Cysteine-rich protein TTG-2) (LIM domain only protein 2) (LMO-2) (T-cell translocation protein 2), LMO2, RBTN2 RBTN1 RHOM2 TTG2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	LMO2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.