



LRRIQ3 Antibody, FITC conjugated

Product Code	CSB-PA013178LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6PVS8
Immunogen	Recombinant Human Leucine-rich repeat and IQ domain-containing protein 3 protein (538-607AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Leucine-rich repeat and IQ domain-containing protein 3 (Leucine-rich repeat-containing protein 44), LRRIQ3, LRRC44
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	LRRIQ3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.