



LYN Antibody, Biotin conjugated

Product Code	CSB-PA013259LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07948
Immunogen	Recombinant Human Tyrosine-protein kinase Lyn protein (135-282AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein kinase Lyn (EC 2.7.10.2) (Lck/Yes-related novel protein tyrosine kinase) (V-yes-1 Yamaguchi sarcoma viral related oncogene homolog) (p53Lyn) (p56Lyn), LYN, JTK8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LYN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.