



MAX Antibody, FITC conjugated

Product Code	CSB-PA013527LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61244
Immunogen	Recombinant Human Protein max protein (1-58AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein max (Class D basic helix-loop-helix protein 4) (bHLHd4) (Myc-associated factor X), MAX, BHLHD4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MAX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.