



ME1 Antibody, HRP conjugated

Product Code	CSB-PA013632LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48163
Immunogen	Recombinant Human NADP-dependent malic enzyme protein (161-519AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADP-dependent malic enzyme (NADP-ME) (EC 1.1.1.40) (Malic enzyme 1), ME1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ME1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.