



ME2 Antibody, Biotin conjugated

Product Code	CSB-PA013633HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23368
Immunogen	Recombinant Human NAD-dependent malic enzyme, mitochondrial protein (296-424AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NAD-dependent malic enzyme, mitochondrial (NAD-ME) (EC 1.1.1.38) (Malic enzyme 2), ME2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ME2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.