



MPL Antibody, HRP conjugated

Product Code	CSB-PA014755LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40238
Immunogen	Recombinant Human Thrombopoietin receptor protein (26-491AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thrombopoietin receptor (TPO-R) (Myeloproliferative leukemia protein) (Proto-oncogene c-Mpl) (CD antigen CD110), MPL, TPOR
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MPL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.