



# MPL Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA014755LC01HU
<b>Abbreviation</b>	Thrombopoietin receptor
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P40238
<b>Immunogen</b>	Recombinant Human Thrombopoietin receptor protein (26-491AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Relevance</b>	Receptor for thrombopoietin. May represent a regulatory molecule specific for TPO-R-dependent immune responses.
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Thrombopoietin receptor (TPO-R) (Myeloproliferative leukemia protein) (Proto-oncogene c-Mpl) (CD antigen CD110), MPL, TPOR
<b>Species</b>	Human
<b>Research Area</b>	Immunology
<b>Target Names</b>	MPL