



MT-CO1 Antibody, Biotin conjugated

Product Code	CSB-PA015072LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00395
Immunogen	Recombinant Human Cytochrome c oxidase subunit 1 protein (474-513AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome c oxidase subunit 1 (EC 1.9.3.1) (Cytochrome c oxidase polypeptide I), MT-CO1, COI COXI MTCO1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MT-CO1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.