



MYO1G Antibody, FITC conjugated

Product Code	CSB-PA015344LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B0I1T2
Immunogen	Recombinant Human Unconventional myosin-Ig protein (209-389AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Unconventional myosin-Ig [Cleaved into: Minor histocompatibility antigen HA-2 (mHag HA-2)], MYO1G, HA2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYO1G
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.