



NDC80 Antibody, HRP conjugated

Product Code	CSB-PA015601LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O14777
Immunogen	Recombinant Human Kinetochores protein NDC80 homolog protein (1-249AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinetochores protein NDC80 homolog (Highly expressed in cancer protein) (Kinetochores protein Hec1) (HsHec1) (Kinetochores-associated protein 2) (Retinoblastoma-associated protein HEC), NDC80, HEC HEC1 KNTC2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NDC80
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.