



PCNA Antibody, HRP conjugated

Product Code	CSB-PA01567B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12004
Immunogen	Recombinant Human Proliferating cell nuclear antigen protein (1-261AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proliferating cell nuclear antigen (PCNA) (Cyclin), PCNA
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PCNA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.