



NEFL Antibody, Biotin conjugated

Product Code	CSB-PA015688LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P07196
Immunogen	Recombinant Human Neurofilament light polypeptide protein (128-271AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurofilament light polypeptide (NF-L) (68 kDa neurofilament protein) (Neurofilament triplet L protein), NEFL, NF68 NFL
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NEFL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.