



NFAT5 Antibody, HRP conjugated

Product Code	CSB-PA015744LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94916
Immunogen	Recombinant Human Nuclear factor of activated T-cells 5 protein (1422-1531AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear factor of activated T-cells 5 (NF-AT5) (T-cell transcription factor NFAT5) (Tonicity-responsive enhancer-binding protein) (TonE-binding protein) (TonEBP), NFAT5, KIAA0827 TONEBP
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NFAT5