



NOXA1 Antibody, HRP conjugated

Product Code	CSB-PA015963LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UR1
Immunogen	Recombinant Human NADPH oxidase activator 1 protein (150-313AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADPH oxidase activator 1 (NOX activator 1) (Antigen NY-CO-31) (NCF2-like protein) (P67phox-like factor) (p51-nox), NOXA1, P51NOX
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NOXA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.