



NR1D1 Antibody, FITC conjugated

Product Code	CSB-PA016043LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20393
Immunogen	Recombinant Human Nuclear receptor subfamily 1 group D member 1 protein (246-558AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor subfamily 1 group D member 1 (Rev-erbA-alpha) (V-erbA-related protein 1) (EAR-1), NR1D1, EAR1 HREV THRAL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NR1D1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.