



NR1H2 Antibody, Biotin conjugated

Product Code	CSB-PA016045HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55055
Immunogen	Recombinant Human Oxysterols receptor LXR-beta protein (1-461AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oxysterols receptor LXR-beta (Liver X receptor beta) (Nuclear receptor NER) (Nuclear receptor subfamily 1 group H member 2) (Ubiquitously-expressed nuclear receptor), NR1H2, LXRB NER UNR
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NR1H2