





NR2F2 Antibody

Product Code	CSB-PA016057LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24468
Immunogen	Recombinant Human COUP transcription factor 2 protein (123-383AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	COUP transcription factor 2 (COUP-TF2) (Apolipoprotein A-I regulatory protein 1) (ARP-1) (COUP transcription factor II) (COUP-TF II) (Nuclear receptor subfamily 2 group F member 2), NR2F2, ARP1 TFCOUP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NR2F2

Image



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA016057LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA016057LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.