



NR2F2 Antibody, FITC conjugated

Product Code	CSB-PA016057LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24468
Immunogen	Recombinant Human COUP transcription factor 2 protein (123-383AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	COUP transcription factor 2 (COUP-TF2) (Apolipoprotein A-I regulatory protein 1) (ARP-1) (COUP transcription factor II) (COUP-TF II) (Nuclear receptor subfamily 2 group F member 2), NR2F2, ARP1 TFCOUP2
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	NR2F2