



ARAF Antibody, Biotin conjugated

Product Code	CSB-PA01679D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10398
Immunogen	Recombinant Human Serine/threonine-protein kinase A-Raf protein (160-310AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine/threonine-protein kinase A-Raf (EC 2.7.11.1) (Proto-oncogene A-Raf) (Proto-oncogene A-Raf-1) (Proto-oncogene Pks), ARAF, ARAF1 PKS PKS2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ARAF
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.