



PCF11 Antibody, FITC conjugated

Product Code	CSB-PA017603LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O94913
Immunogen	Recombinant Human Pre-mRNA cleavage complex 2 protein Pcf11 protein (409-604AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-mRNA cleavage complex 2 protein Pcf11 (Pre-mRNA cleavage complex II protein Pcf11), PCF11, KIAA0824
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PCF11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.