



PEBP1 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA017766HB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P30086 |
| Immunogen | Recombinant Human Phosphatidylethanolamine-binding protein 1 protein (97-128AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Phosphatidylethanolamine-binding protein 1 (PEBP-1) (HCNPpp) (Neuropolypeptide h3) (Prostatic-binding protein) (Raf kinase inhibitor protein) (RKIP) [Cleaved into: Hippocampal cholinergic neurostimulating peptide (HCNP)], PEBP1, PBP PEBP |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | PEBP1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |