



PEBP1 Antibody, Biotin conjugated

Product Code	CSB-PA017766HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30086
Immunogen	Recombinant Human Phosphatidylethanolamine-binding protein 1 protein (97-128AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylethanolamine-binding protein 1 (PEBP-1) (HCNPpp) (Neuropolypeptide h3) (Prostatic-binding protein) (Raf kinase inhibitor protein) (RKIP) [Cleaved into: Hippocampal cholinergic neurostimulating peptide (HCNP)], PEBP1, PBP PEBP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PEBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.