



PGP Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA017873LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | A6NDG6 |
| Immunogen | Recombinant Human Glycerol-3-phosphate phosphatase protein (140-300AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Glycerol-3-phosphate phosphatase (G3PP) (EC 3.1.3.21) (Aspartate-based ubiquitous Mg(2+)-dependent phosphatase) (AUM) (EC 3.1.3.48) (Phosphoglycolate phosphatase) (PGP), PGP |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | PGP |