



PGP Antibody, Biotin conjugated

Product Code	CSB-PA017873LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A6NDG6
Immunogen	Recombinant Human Glycerol-3-phosphate phosphatase protein (140-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycerol-3-phosphate phosphatase (G3PP) (EC 3.1.3.21) (Aspartate-based ubiquitous Mg(2+)-dependent phosphatase) (AUM) (EC 3.1.3.48) (Phosphoglycolate phosphatase) (PGP), PGP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PGP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.